

## 193nm Narrow Linewidth Fiber Lasers

Based on single frequency seed laser, low noise fiber amplifier and high efficiency frequency convert technique, Precilasers can provide narrow linewidth 193nm CW lasers, features with high frequency stability, long life and good beam quality.

### Features

- Narrow Linewidth
- Good Beam Quality
- Never Mode Hop
- CW Laser

### Application

- Precision Measurement



Specification	
Partnumber	FL-SF-193-X <sup>(1)</sup> -CW
Center Wavelength	193nm
Linewidth <sup>(2)</sup> (100 us integration )	<100kHz
Operational Mode	CW
Output Power	>3mW
Beam Diameter	1±0.1 mm
Beam Waist Location	output aperture±1m
Beam Divergence (full angle)	<0.4mrad
Beam Circularity	>90%
Beam Pointing Stability	<50urad @30min, RMS
Beam Quality	TEM <sub>00</sub> , M <sup>2</sup> <1.5
Polarization	linear, >20dB
Power Stability	<1% @2hours, RMS
Relative Intensity Noise	<0.5%(RMS, 10Hz-10MHz integration)
Cooling	Water cooling

Warm-Up Time	
Cold Start	<2 hours
Warm Start	<30 minutes

Demension and Weight	
Controller Size	≤483mm x 420mm x 133mm ( 3mW )
Laser Head Size	≤550mm x 350mm x 131mm ( 3mW )
Chiller Size	≤58cm x 29cm x 52cm
Laser Head Weight	≤30kg
Controller Weight	≤30kg
Chiller Weight	≤25kg
Cable Length Between Laser Head And Controller	>10m

Environment and Power Supply	
Operating Temperature ( No condensing )	15-35 °C
Power Supply	100V-240V, AC, 50/60Hz
Power Consumption	<400W
Communication	RS422


Option	
Frequency Stabilization Option	Wavelength stability<±10fm@24h
Purge Air Option	Provide air pump

(1)X indicates laser power,unit in W.For example output power is 3mW : FL-SF-193-0.003-CW.


(2)Linewidth measured by delayed self-heterodyne

**Water cooling version size  
(Output power: 3mW)**

Laser Controller



laser Head



Non-dismantle fiber length  $\geq 10\text{m}$  between controller and laser head



Shanghai Precilasers Technology Co., Ltd.  
 Floor 2, Building 2, No. 1918, Xupan Road, Jiading  
 District, Shanghai  
 021-59160265

[info@Precilasers.com](mailto:info@Precilasers.com) | [www.Precilasers.com](http://www.Precilasers.com)



**△ Laser Hazard**

Visible or invisible laser radiation, avoid eye or skin exposure to direct, reflected or filtered radiation.

CLASS 4 Laser Products

